Instructions for use

prodisc® L – Lumbar disc prosthesis

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pro**disc L** Lumbar disc prosthesis

Safety precautions

Please read these instructions for use and the corresponding surgical techniques carefully before use. Ensure that you are familiar with the appropriate surgical technique.

Contents

The pro**disc L** Lumbar disc prosthesis is made up of three components:

- pro**disc L** superior endplate
- pro**disc L** inferior endplate
- pro**disc L** inlay

All implant components are packaged together using a double sterile barrier method.

Description

The components of the pro**disc L** Lumbar disc prosthesis are made from:

- (1) Superior and inferior endplate: CoCrMo (Co-28Cr-6Mo) per ISO 5832-12 with pure titanium coating per ISO 5832-2
- (2) Inlay: UHMWPE per ISO 5834-2 with an integrated tantalum x-ray marker per ISO 13782.

Contents are supplied sterile.

English

Intended use

prodisc L implants are used to replace a lumbar intervertebral disc and to restore disc height and segmental motion.

Indications

Lumbar discopathy at one or two contiguous levels.

Contraindications

- Spinal stenosis, radiculopathy
- Increased segmental instability
- Spinal deformities, spondylolisthesis above 25%
- Radiographic confirmation of severe facet joint disease or degeneration
- Osteoporosis, osteochondrosis, and severe osteopenia
- Acute or chronic systemic, spinal, or localized infections
- Systemic and metabolic diseases
- Any medical and surgical conditions precluding the potential benefit of spinal surgery
- Foreign body sensitivity to the implant materials
- Dependency on pharmaceutical drugs, drug abuse, or alcoholism
- Pregnancy
- Obesity
- Lack of patient cooperation
- Skeletally immature patients

Patient exclusion recommendations

Patient selection is one of the most important factors contributing to the outcome of the total disc replacement procedure. The following may affect clinical outcomes:

- A condition of senility or mental illness, alcoholism or smoking
- Dependency on pharmaceutical drugs or drug abuse
- The patient's occupation or activity level
- Compromised vertebral bodies at affected level due to current or past trauma (fractures)
- Substantial loss of disc height, where applied segmental distraction may lead to damage of the great vessels
- Involved vertebral endplate dimensionally smaller than the minimum implant footprint size in both the medial-lateral and the anterior-posterior directions
- Severe abnormality of the endplate (e.g. large Schmorl nodes)
- Posterior spinal defect (e.g. Pars defect)

Precautions

- Proper surgical performance of the implantation is the responsibility of the operating surgeon.
- The operating surgeon must have a thorough command of both the hands-on and conceptual aspects of the established operating techniques.
- Assembling and implanting the implant components is the responsibility of the operating surgeon.
- The manufacturer is not responsible for any complications arising

English

from incorrect diagnosis, choice of incorrect implant, incorrectly combined implant components and/or operating techniques, the limitations of treatment methods, or inadequate asepsis.

- Under no circumstances may modular implant components from different suppliers be combined.
- The implant components applied (name, article number, lot number) must be documented in each patient's record.
- During the postoperative phase, in addition to mobility and muscle training, it is of particular importance that the physician keeps the patient well informed.
- Damage to the weight-bearing structures can give rise to loosening of the components, dislocation and migration, as well as to other grave complications. To ensure the earliest possible detection of such catalysts of implant dysfunction, the lumbar disc prosthesis must be checked periodically post operative, using appropriate techniques.

Potential risks

Potential risks to health associated with the use of the pro**disc L** include, but are not limited to:

- Those commonly associated with any surgical procedure
- Death as a potential risk associated with any surgical procedure
- Those specifically associated with lumbar spinal surgery using an anterior approach
- Those associated with a lumbar disc prosthesis (including the pro**disc L**)
- Bending, loosening, wear and tear on, or fracture of implant components

- Loss of fixation, dislocation and migration
- Neurological injury
- Injury to vessels, nerves and organs
- Primary and secondary infections
- Allergic reactions to implant materials
- Tissue reaction to implant materials
- Venous thrombosis, lung embolism and cardiac arrest
- Hematoma and impaired wound healing
- Periarticular calcification and fusion
- Subsidence
- Vertebral body fracture
- Spine alignment changes
- Degenerative changes of the index or adjacent level

Procedure

- The **prodisc L** Lumbar disc prosthesis must be implanted with the **prodisc L** trial implant with adjustable stop and **prodisc L** instruments only.
- The **prodisc L** trial implant with adjustable stop and **prodisc L** instruments can be ordered from the manufacturer at any time.

The operating surgeon draws up an operation plan specifying and documenting the following:

- Implant components and their dimensions.
- Proper position of the implant components in the intervertebral space.
- Determination of intraoperative orientation points.

English

The following conditions must be fulfilled prior to application:

- All requisite implant components are readily available.
- Highly aseptic operating conditions are present.
- All requisite implantation instruments must be available and in working order, including **prodisc L** trial implant with adjustable stop and **prodisc L** instruments.
- Damaged or defective instruments should not be used or processed. Contact your local Centinel Spine representative or dealer for repair or replacement.
- The use of an instrument for tasks other than those for which they are intended may result in damaged / broken instruments or patient injury.
- The operating surgeon and operating room team must be thoroughly familiar with the operating technique, as well as the range of implants and instruments to be applied; complete information on these subjects must be readily available at the workplace.
- The operating surgeon must be especially trained in spinal surgery, biomechanical principles of the spine and the relevant operating techniques.

The operative procedure has been explained to the patient, and the patient's understanding of the following information has been documented:

- The patient is aware of the risks associated with neurosurgery, general surgery, orthopaedic surgery and with general anaesthesia.

- The patient has been informed about the advantages and disadvantages of a lumbar disc prosthesis and about possible alternative treatments.
- The lumbar disc prosthesis can fail owing to excessive load, wear and tear, or infection.
- The service life of the implant is determined by body weight.
- The lumbar disc prosthesis must not be subjected to overload through extreme strain, or through work-related or athletic activities.
- Corrective surgery may be necessary if the implant fails.
- In the event corrective surgery is performed, it may not be possible to restore segmental motion.
- The patient must have their physician carry out follow-up examinations of the lumbar disc prosthesis at regular intervals.

Disinfection, cleaning, care and sterilization

- The implant components come individually packed in protective packaging that is labeled according to its contents.
- The implant components are radiosterilized (gamma sterilization, minimum 25 kGy).
- Store implant components in their original protective packaging, and do not remove them from the packaging until immediately before use.
- Prior to use, check the product expiration date and verify the integrity of the sterile packaging.
- Do not use implant components that are past their expiration dates or whose packaging is damaged.
- Never use either damaged implant components or implant compo-

English

nents that have been removed surgically.

- Never resterilize or reuse implants that have a coated surface if such implants have been contaminated by blood or secretions.
- Never resterilize implants either made of polyethylene, or implants containing polyethylene components.
- The implants are not designed to be disinfected or cleaned by the user.

Magnetic Resonance environment

MR Conditional:

Non-clinical testing of the worst-case scenario has demonstrated that the implants of the **prodisc L** system are MR conditional. These articles can be scanned safely under the following conditions:

- Static magnetic field of 1.5 Tesla and 3.0 Tesla
- Spatial gradient magnetic field of 90 mT/cm (900 Gauss/cm)
- Maximum whole body averaged specific absorption rate (SAR) of 2 W/kg for 15 minutes of scanning

In non-clinical testing, the **prodisc L** implant will produce a temperature rise not greater than 1.8 °C at a maximum whole body averaged specific absorption rate (SAR) of 2 W/kg, as assessed by calorimetry for 15 minutes of MR scanning in a 1.5 Tesla and 3.0 Tesla MR scanner.

MR Imaging quality may be compromised if the area of interest is

in the exact same area or relatively close to the position of the **prodisc L** device.

Warranty

All warranty rights are lost if repairs or modifications are carried out by an unauthorized service center. The manufacturer does not take responsibility for any effects on safety, reliability or performance of the product if the product is not used in conformity with the instructions for use.

For further information

If further information on this product is needed, please contact your local Centinel Spine representative or dealer.

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trademark of Centinel Spine, LLC.
affiliates. All rights reserved.
CENTINEL SPINE GmbH
Centinel Spine Schweiz GmbH
Centinel Spine Australia Pty Ltd
© 2021 Centinel Spine, LLC or its