Instructions for Use

Prodisc-Oblique – Lumbar Disc Prosthesis
using the Antero-Lateral Approach

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Prodisc-Oblique – lumbar disc prosthesis using the antero-lateral approach

Safety precautions

Please read these instructions for use, the Synthes brochure “Important Information” and the corresponding surgical techniques carefully before use. Ensure that you are familiar with the appropriate surgical technique.

Contents

The three components of the Prodisc-Oblique lumbar disc prosthesis are supplied separately and consist of the following items:

- Prodisc-Oblique superior endplate
- Prodisc-Oblique inferior endplate
- Prodisc-Oblique PE-inlay

Description

The three components of the Prodisc-Oblique lumbar disc prosthesis are made of the following materials:

- (1) Superior and inferior endplate: CoCrMo (Co-28Cr-6Mo) per ISO 5832-12 with pure titanium coating per ISO 5832-2.
- (2) UHMWPE Inlay per ISO 5834-2 with an integrated tantalum x-ray marker per ISO 13782.

Contents are supplied sterile.

Intended use

Prodisc-Oblique implants are used to replace lumbar intervertebral discs and to restore disc height and segmental motion at vertebral levels L1 to L5. They are inserted by using the antero-lateral approach.

Indications

Lumbar discopathy for levels L1 to L5, for which patient's anatomy allows the antero-lateral approach.

Contraindications

- Vertebral level L5/S1
- Any circumstances preventing antero-lateral access to the spinal column (including extensive abdominal and/or retroperitoneal surgery from the left)
- Compromised vertebral bodies at affected level due to current or past trauma (fractures)
- Substantial loss of disc height, where applied segmental distraction may lead to damage of the great vessels
- Isolated radicular compression syndromes, especially due to disc herniation
- Bony lumbar spinal stenosis
- Spondylolysis/Retroisthesis
- Pars defect
- Complete laminectomy
- Osteoporosis and osteopenia
- Active systemic infection or infection localized to the site of implantation

- Allergy or sensitivity to implant materials (foreign body sensitivity to the implant materials)
- Adipositas
- Pregnancy
- Involved vertebral endplate dimensionally smaller than the minimum implant footprint size in both the medial-lateral and the anterior-posterior directions
- Severe abnormality of the endplate (e.g. large Schmorl nodes)
- Facet joint disease or degeneration
- Systemic and / or metabolic diseases
- Back or leg pain of unknown etiology
- Active malignancy (e.g. tumors)
- Acute or chronic infections (systemic and / or local)
- Dependency on pharmaceutical drugs or drug abuse, or alcoholism
- Predominant psychosocial factors / illnesses
- Foraminal or lateral spinal canal stenosis
- Lack of patient compliance
- Any cases not listed in the indications

Patient exclusion recommendations

Patient selection is one of the most important factors contributing to the outcome of the total disc replacement procedure. The following may affect clinical outcomes:

- A condition of senility or mental illness, alcoholism or smoking
- Dependency on pharmaceutical drugs or drug abuse
- The patient's occupation or activity level

English

- Compromised vertebral bodies at affected level due to current or past trauma (fractures)
- Substantial loss of disc height, where applied segmental distraction may lead to damage of the great vessels
- Involved vertebral endplate dimensionally smaller than the minimum implant footprint size in both the medial-lateral and the anterior-posterior directions
- Severe abnormality of the endplate (e.g. large Schmorl nodes)

Precautions

- Proper surgical performance of the implantation is the responsibility of the operating surgeon.
- The operating surgeon must have a thorough command of both the hands-on and conceptual aspects of the established operating techniques.
- Assembling and implanting the implant components is the responsibility of the operating surgeon.
- The manufacturer is not responsible for any complications arising from incorrect diagnosis, incorrect choice of implant, incorrectly combined implant components and/or operating techniques, the limitations of treatment methods, or inadequate asepsis.
- Modular implant components from different manufacturers may not be combined.
- The implant components applied (name, reference number, lot number) must be documented in the patient record.
- During the postoperative phase, in addition to mobility and muscle training, it is of particular importance that the physician keep the patient well informed.

- It is mandatory that the patient actively cooperate and follow the instructions given by the physician.
- Damage to weight-bearing structures can give rise to loosening of the components, migration, as well as to other grave complications. To ensure the earliest possible detection of such catalysts of implant dysfunction, the lumbar disc prosthesis must be checked periodically post operatively, using appropriate techniques.

Potential risks

Potential risks to health associated with the use of the Prodisc-Oblique include, but are not limited to:

- Those commonly associated with any surgical procedure
- Death as a potential risk associated with any surgical procedure
- Those specifically associated with lumbar spinal surgery using an antero-lateral approach
- Those associated with a lumbar disc prosthesis (including the Prodisc-Oblique)
- Bending, loosening, wear and tear on, or fracture of the implant components
- Dislocation, loss of fixation and migration
- Neurological injuries
- Injury to vessels, nerves and organs
- Primary and secondary infections
- Allergic reactions to the implant materials
- Tissue reactions to the implant materials
- Venous thrombosis, pulmonary embolism and cardiac arrest
- Haematoma and impaired wound healing

English

- Periarticular calcification and fusion
- Subsidence
- Vertebral body fracture
- Spine alignment changes
- Degenerative changes of the index or adjacent level

Procedure

The Prodisc-Oblique may only be implanted with the instruments provided for this purpose.

The surgeon draws up an operative plan, specifying and recording the following:

- Implant components and their dimensions
- Proper positioning of the implant components within the bone
- Determination of the intraoperative orientation points

The following conditions must be fulfilled prior to application:

- All requisite implant components must be readily available.
- Highly aseptic operating conditions are present.
- All requisite instruments for the implantation must be sterile, with all components available and in good working order.
- Damaged or faulty instruments may not be used or processed. Contact your local Synthes representative or dealer for repair or replacement.
- The use of an instrument for tasks other than those for which they are intended may result in damaged / broken instruments.

- The operating surgeon and operating room team must be familiar with the operating technique, as well as the range of implants and instruments to be applied.
- The operating surgeon must be especially trained in spinal surgery, bio-mechanical principles of the spine and the relevant operating techniques.

The operative procedure has been explained to the patient, and the patient's understanding of the following information has been documented:

- The patient is aware of the risks associated with neurosurgery, general surgery, orthopedic surgery and with general anesthesia.
- The patient has been informed of the advantages and disadvantages of a lumbar disc prosthesis and possible alternative treatments.
- The Prodisc-Oblique disc prosthesis can fail owing to excessive load, wear and tear or infection.
- The expected useful life of the implant is determined by body weight and activity.
- The Prodisc-Oblique lumbar disc prosthesis must not be subjected to overload through extreme strain, or through work-related or athletic activities.
- Corrective surgery may be necessary if the implant fails.
- In the event that corrective surgery is performed it may not be possible to restore segmental motion.
- The patient must have their physician carry out follow-up examinations of the lumbar disc prosthesis at regular intervals.
- The patient must follow the postoperative treatment methods that the surgeon stipulates.

Disinfection, cleaning, care and sterilization

- The implant components are supplied individually packed in a protective packaging that is labeled according to its contents.
- The implant components are radiosterilized (gamma sterilization, minimum 25 kGy).
- Store implant components in their original protective packaging, and do not remove them from the packaging until immediately before use.
- Prior to use, check the product expiration date and verify the integrity of the sterile packaging.
- Do not use implant components that are past their expiration dates or whose packaging is damaged.
- Never use either damaged implant components or implant components that have been removed surgically.
- Never resterilize or reuse implants that have a coated surface if such implants have been contaminated by blood or secretions.
- Never resterilize implants either made of polyethylene, or implants containing polyethylene components.
- Never resterilize or reuse the Prodisc-Oblique lumbar disc prosthesis. The implants are not designed to be cleaned or disinfected by the user.

Magnetic Resonance environment

MR Conditional:

Non-clinical testing of the worst-case scenario has demonstrated that the implants of the **Prodisc-Oblique** system are MR conditional. These articles can be scanned safely under the following conditions:

- Static magnetic field of 1.5 Tesla and 3.0 Tesla.
- Spatial gradient field of 90 mT/cm (900 Gauss/cm).
- Maximum whole body averaged specific absorption rate (SAR) of 2 W/kg for 15 minutes of scanning.

Based on non-clinical testing, the **Prodisc-Oblique** implant will produce a temperature rise not greater than 1.8 °C at a maximum whole body averaged specific absorption rate (SAR) of 2 W/kg, as assessed by calorimetry for 15 minutes of MR scanning in a 1.5 Tesla and 3.0 Tesla MR scanner.

MR Imaging quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the **Prodisc-Oblique** device.

Warranty

All warranty rights are invalidated if repairs or modifications are carried out by an unauthorized service center. The manufacturer accepts no liability for any impact on safety, reliability or performance of the product if the product is not used in conformity with these instructions for use.

For further information

If further information on this product is needed, please contact your local Synthes representative or dealer.

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© DePuy Synthes Spine, a division of Synthes GmbH. 2016. All rights reserved.
Synthes GmbH